

IMMUNO-DOT BLOT ASSAY

(Wright & Morton, 1989)

INTRODUCTION

This assay may be used as a positive or negative test for presence of glomalin but not for determining concentration. It also gives a permanent result, since the color is on the nitrocellulose membrane and does not fade.

MATERIALS

0.2 um nitrocellulose membrane
petri dish or dish to incubate the membrane in
2% non-fat milk (2g powdered milk/100ml PBS)
PBS (phosphate buffered saline)
PBST (PBS with Tween 20)
1% BSA (1 g BSA/100 mls PBS) (Make about 500-1000 mls of stock and dispense in 6 ml aliquots that will be frozen until needed.)
TBS (Tris buffered saline (20mM tris and 250mM NaCl)
MAb32B11 antibody (diluted with PBS, typically 1 ml in 5 mls PBS)
Biotinylated anti-mouse IgM (4.8 ul/6 ml 1% BSA)
ExtrAvidin peroxidase (3.0 ul/ 6 ml 1% BSA)
Color developer (make just prior to using)
 0.015g 4-chlor-1-naphthol in 5ml ice cold methanol,
 right before use add 25 ml TBS + 15ul of 30% hydrogen peroxide
dissecting needle
shaker or tilt table

METHODS

- (1) Mark the membrane into small squares with the dissecting needle
- (2) Place 1ul of undiluted sample in each square
- (3) Block membrane with 2% non-fat milk by shaking for 15 min.
- (4) Add diluted MAb32B11 and incubate on shaker for 1 hr.
- (5) Remove antibody, add PBST and incubate on shaker for 5 min. Repeat PBST incubation twice.
- (6) Incubate diluted biotinylated anti-mouse IgM with the membrane for 1 hr on shaker.
- (7) Remove IgM, add PBST and incubate on shaker for 5 min. Repeat PBST incubation twice.
- (8) Add diluted ExtrAvidin peroxidase and incubated for 1 hr on shaker.
- (9) Remove peroxidase, add PBST and incubate on shaker for 5 min. Repeat PBST incubation twice and once with TBS.
- (10) Develop with color developer until color is seen. Remove and dry. Store dry at room temp.